



POWER SYSTEMS

AESCULAP® Acculan 4 Electro

LOOK. FEEL. EXPERIENCE.

AESCULAP® Acculan 4 Electro

LOOK. FEEL. EXPERIENCE.



Acc

Acculan 4 Electro offers an alternative to battery-driven power systems for the use in Orthopaedics, Traumatology, Cardiothoracic Surgery and Dermatology without compromises in power or reliability.

THE DIFFERENCE IS IN THE DETAIL

AESCULAP® pioneered control electronics that are not exposed to the stress of reprocessing. Our Electro solution, available for all Acculan 4 handpieces, continues this proven concept.

MECHANICAL PROCESSING

Thanks to the titanium housing and the ECCOS storage solution Acculan 4 is suitable for mechanical cleaning, even with alkaline cleaning media. With our rinsing adapter, we offer you a simple way for cleaning the inside of lumen mechanically.

SERVICE LABEL

Regular servicing and preventive maintenance is essential to ensure the reliability of motor systems over a long service life. The service label on each Acculan 4 machine clearly shows when the next preventive maintenance is due.

Acculan 4

Combine the Acculan 4 handpieces ...

GA330 Acculan 4 drill and reamer

GA331 Acculan 4 oscillating saw

GA334 Acculan 4 reciprocating saw

GA340 Acculan 4 dermatom 0.2 mm

GA344 Acculan 4 small drill


... with the Electro components

GA875 Acculan 3Ti Electro lid

GA876 Acculan Electro adapter

GA877 Acculan 3Ti Electro console

GA878 Acculan 3Ti Electro connection cable



- Accessories, such as the ECCOS storage solution, can be found in the Power Systems main catalogue 084002.
- The tools (e.g. dermatome blades or skin carriers) can be found in the Burrs & Blades catalogue 017599.

AESCULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" is a registered trademark of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.